

CAPPING - DECAPPING

ErgoCap Tube Decapper / Recapper

The ErgoCap provides an easy and pain-free method to uncap and cap a variety of vials and tubes. To use, simply place the vial into the opening with minimal pressure and the ErgoCap will do the rest, either uncapping or recapping the vial. The ErgoCap will be indispensable for routine manual operations that require constant decapping and recapping and will minimize the risk of repetitive motion injuries in the lab.



- Works with a variety of tube sizes
- Dimensions 110 mm W x 338mm D x 295 mm H
- Weight 5.5 kgs (12 lbs)
- 100 -0240 VAC power compatible

