



SEALING

Heat Seal Rolls for Automated Heat Sealers



Biosero offers an extensive range of sealing instruments and consumables, including heat seals compatible with the major thermal sealers on the market, including seals specifically designed for use with a wide variety of solvents including DMSO and "crystal-clear" films for protein crystallography and QPCR applications.

Heat seals are compatible with models from KBiosystems, Abgene (Thermo) and Velocity 11 and are available at competitive prices.



	Peelable Foil	Pierceable Foil	Clear Peelable	Heavy Duty Foil
ALPS 300™ / ALPS 3000™ & Wasp	500127 (610m)	500131 (610m)	500119 (500m)	500141 (610m)
	500128 (25m)	500132 (25m)	500136 (25m)	500142 (25m)
			500139 (610m)	500145 (500m)
			500140 (25m)	500146 (25m)
			500143 (610m)	
PlateLoc™	500129 (500m)	500133 (500m)	500137 (500m)	500147 (500m)
	500130 (25m)	500134 (25m)	500138 (25m)	500148 (25m)
	Clear Weld Seal (cannot be peeled, difficult to pierce)	DMSO Resistant Peelable	Clear Pierceable for ABI 3730	
PlateLoc™	500159 (500m)	500162 (500m)	500167 (500m)	
		5500163 (25m)	500138 (25m)	

If you require a seal type that is not shown here, contact us to check for availability.

