

SEALING

KAPS-500 Auto Adhesive Plate Sealer



- Desktop Automated Cold/Adhesive Plate Sealing System
- Seal SBS format micro well plates
- Seal plate heights from 9 to 48mm high
- Variable sealing pressure
- Variable sealing speed
- Easy change pressure rollers for different film types requirements
- Plate holders supplied for different plate types support requirements due to top surface pressure
- Seal time from 4 plates per minute
- Foil low sensing
- Touch panel display
- Screen or foot pedal manual operation
- System only 25kgs (52lbs)
- Stacker Compatible with KBiosystems Scorpion stacker for increased throughput
- OEM Stacker compatible call for information
- Integrator Friendly
- 110/220 vac power compatible